


EEMB DS Capstone Tracking, Key Values, and Timestamping GUI Manual







## **EEMB DS Capstone** Home Obtain Key Stats + Timestamp Tracking Video: /Users/ashleyson/Desktop/tria Select Video Background: /Users/ashleyson/Desktop/dra Select Background Output File Path: /Users/ashleyson/Desktop/dra Select Output Path Start Second: 30 Run Tracking

Once you have selected a video background, output file path, and start second, click "Run Tracking".



Once you have ran the tracking, there will be a pop up for you to initialize the four corners of the tank. Select the four corners. If you are not content with the four corners, click "n" on your keyboard to clear the points for you to reselect. If you are satisfied with the four corners, click "y" and it will direct you to the next steps.









Once selecting "OK," the tracking process will begin, you can find the output.csv in the same location as your video and background.

Happy Tracking! ☺






Input the closeness threshold. This value is measured in pixels. Then click on the "Input Closeness Threshold:" button.



Input the identification name for the 3 csv files. Then click on the "Input File Identification Name:" button.



View above

Input the number of objects to track. Then click on the "Input Amount of Objects to Track:" button.

Once selecting "Generate Important Timestamps from this Dataset," the process will begin, you can find the csv files in your designated output path.